

# ExFuse™

Demineralized Bone Matrix with Cancellous Bone



DFDBA



Cancellous



Gel



Putty



DFDBA



Cancellous



Cat No.	Size	Particle Size	Type		
EG03	0.3cc	200~850µm	Gel		
EG05	0.5cc				
EG1	1cc				
EG3	3cc				
EG5	5cc				
EG8	8cc				
EG10	10cc				
EP03	0.3cc			200~850µm	Putty
EP05	0.5cc				
EP1	1cc				
EP3	3cc				
EP5	5cc				
EP8	8cc				
EP10	10cc				



DFDBA



Cancellous

# ExFuse™ Gel/Putty

## 1. COMPOSITION

- 1-1. Main components: Demineralized Cortical Bone
- 1-2. Additives: Cancellous Bone
- 1-3. Additives: CMC (Carboxymethyl Cellulose Carrier)  
CMC is used as additive intends to make DBM easier to handle.

## 2. PRODUCT TYPE

- 2-1. DBM Gel  
DBM Gel is thinner than Putty type and can be directly applied into the defect site since it is provided in a syringe.
- 2-2. DBM Putty  
DBM Putty is available in a syringe and vial and has a high handle-ability for ease of use.

## 3. SIZE

- |                                   |                                   |
|-----------------------------------|-----------------------------------|
| 3-1. DBM Gel                      | 3-2. DBM Putty                    |
| 0.3 / 0.5 / 1 / 3 / 5 / 8 / 10 cc | 0.3 / 0.5 / 1 / 3 / 5 / 8 / 10 cc |



## 4. INDICATION

- 4-1. Orthopedic / Spine surgery
  - Compensation and reconstruction of autograft bone.
  - Reconstruction of maxillofacial defect
  - Spine fixation with adequate materials
  - Bone fixation and filling in various sizes and forms of bone defect
  - Spinal fusions
  - Oncology
  - Joint revision
  - Long bone trauma
  - Distal fractures
  - Small bone procedure
  - Craniomaxillofacial
- 4-2. Dental surgery
 

<p><b>Implantological Indications</b></p> <ul style="list-style-type: none"> <li>- Extraction Site</li> <li>- Ridge Preservation</li> <li>- Sinus Floor Elevation</li> <li>- Socket Preservation</li> <li>- Wall Defect</li> <li>- Torus</li> <li>- Fenestration Defect</li> <li>- Dehiscence Defect</li> </ul>	<p><b>Periodontal Indications</b></p> <ul style="list-style-type: none"> <li>- Periodontal Augmentation</li> <li>- Exposed Implant</li> <li>- Periodontal Defect</li> <li>- Esthetic Augmentation</li> <li>- Cysts Tumors</li> </ul>
---	--

## 6. PACKAGE

### 6-1. Double-sealed Blister




### 6-2. Outer Paper Package




### 6-3. Insert

- Instruction for Use
- Transplantation Record
- Inner Label

### 6-4. How to Supplied



Gel Type



Putty Type